

A.4513

A.4513I

A.4520

A.4520I

A.4525

A.4525I

A.4525V

A.4525VI

Rectoscopy kits. Sterile.

sapi med

DIAGNOSTIC
LINE



CE
0476

AVAILABLE IN NON STERILE VERSIONS AS WELL.

TECHNICAL FEATURES

- Flash free
- Guaranteed atraumatic
- Compatible with all Sapi Med light sources
- External cm markings on tubes
- They are made in plastic material conforms to the legislation 10993-1 regarding the biocompatibility
- Disposable grip
- Latex free
- Sterile/Non sterile

HANDGRIP

Order code	Description	Quantity per box
	REF A.4513 Black polymer handgrip with 130mm proctoscope. Disposable & Sterile .	25
	REF A.4520 Black polymer handgrip with 200mm rectoscope. Disposable & Sterile .	25
	REF A.4525 Black polymer handgrip with 250mm rectoscope. Disposable & Sterile .	25
	REF A.4525V Black polymer handgrip with 250mm green rectoscope. Disposable & Sterile .	25
	REF A.4513NS Black polymer handgrip with 130mm proctoscope. Disposable & Non Sterile .	25
	REF A.4520NS Black polymer handgrip with 200mm rectoscope. Disposable & Non Sterile .	25
	REF A.4525NS Black polymer handgrip with 250mm rectoscope. Disposable & Non Sterile .	25
	REF A.4525VNS Black polymer handgrip with 250mm green rectoscope. Disposable & Non Sterile .	25

HANDGRIP WITH INSUFFLATOR

Order code	Description	Quantity per box
	REF A.4513I Black polymer handgrip with a single use insufflator & 130mm proctoscope. Disposable & Sterile .	20
	REF A.4520I Black polymer handgrip with a single use insufflator & 200mm rectoscope. Disposable & Sterile .	20
	REF A.4525I Black polymer handgrip with a single use insufflator & 250mm rectoscope. Disposable & Sterile .	20
	REF A.4525VI Black polymer handgrip with a single use insufflator & 250mm green rectoscope. Disposable & Sterile .	20
	REF A.4513INS Black polymer handgrip with a single use insufflator & 130mm proctoscope. Disposable & Non Sterile .	20
	REF A.4520INS Black polymer handgrip with a single use insufflator & 200mm rectoscope. Disposable & Non Sterile .	20
	REF A.4525INS Black polymer handgrip with a single use insufflator & 250mm rectoscope. Disposable & Non Sterile .	20
	REF A.4525VINS Black polymer handgrip with a single use insufflator & 250mm green rectoscope. Disposable & Non Sterile .	20