

A.4113

*Disposable medical devices
for rectosigmoidoscopy.*

sapi med
DIAGNOSTIC
LINE

A.4120

A.4125V

A.4125

A.4131

A.8200




A complete range of single-use, diagnostic and paediatric rectoscopes with multi-use or singleuse grip equipped with a magnifying glass and an attachment for air insufflator. Easy to use, handy, perfectly smooth, atraumatic.

TECHNICAL FEATURES

- Flash-free
- Guaranteed atraumatic
- Compatible with all Sapi Med light sources
- External cm markings on tubes
- They are made in plastic material conforms to the legislation 10993-1 regarding the biocompatibility
- Latex free
- Non sterile

REF A.4113 €€



Short rectoscope

Transparent tube length: 130mm


Internal diameter: 17mm

External diameter: 19,6mm

Quantity per box: 25 scopes

Kit (25 scopes + 2 grips): REF A.4300

REF A.4120 €€



Paediatric rectoscope

Transparent tube length: 200mm


Internal diameter: 12,5mm

External diameter: 15,1mm

Quantity per box: 25 scopes

Kit (25 scopes + 2 grips): REF A.4400

REF A.4125 / A.4125V €€



Rectoscope

A.4125: Transparent tube

A.4125V: Green tube

Length: 250mm

Internal diameter: 17mm

External diameter: 19,6mm

Quantity per box: 25 scopes

Kit (25 trasp. scopes + 2 grips): REF A.4200

Kit (25 green scopes + 2 grips): REF A.4200V

REF A.4130 €€

Still available rigid-flexible rectosigmoidoscope 300mm long.

REF A.4131 €€



Grip


Tough, multi-use plastic grip. For models:

A.4113 - A.4120 - A.4125 - A.4125V - A.4130

Equipped with:

Hinged closure incorporating lens and insufflator connection. Easy-toattach press fitting for cold light source cable. Cleaning instructions for limited re-use.

REF A.8200 €€ / A.8200ST (sterile) €€



Disposable insufflator bellow for air insufflations.

Material: in soft PVC

Insufflating capacity: 65cm³

Weight: 60g

Latex free

Quantity per box: 10 pcs. or 50 pcs.