

A.4018L

Self·led·light®

sapi med
DIAGNOSTIC
LINE

A.4019L

A.4021L

A.4022L

A.4023L

A.4024L



CE

Disposable **Self-led-light** anoscopes can be used by the specialist both in clinic and in the operating theatre.

Made of transparent plastic material and supplied in individually sealed bags with led light built-in. Led lighting can be performed by pulling out the transparent strip. Guaranteed same light intensity for two hours usage.


Compatible, through adapter, to systems of image aquisition (cameras) **ALREADY AVAILABLE IN HOSPITALS** (for example operating theatres laparoscopic column).

ALL ANOSCOPES WITH BUILT-IN LED LIGHT ARE AVAILABLE IN THE STERILE VERSION.

TECHNICAL FEATURES


- Clear, transparent plastic proctoscopes
- Flash free
- Guaranteed atraumatic
- Not rechargeable led lighting
- Separable anoscope and led lighting for safe disposable
- Conforms fully with essential requirements of European Council Directive 93/42/EU
- Latex free
- Non sterile/Sterile

REF A.4018L CE / A.4018LST (sterile) CE




Proctoscope paediatric/diagnostic flute beaked with led lighting built-in
 Tube length: 43/56mm
 Tube diameter: 18mm
Quantity per box: 25 pcs.

REF A.4019L CE / A.4019LST (sterile) CE




Proctoscope paediatric/diagnostic trunk-end with led lighting built-in
 Tube length: 54mm
 Tube diameter: 18mm
Quantity per box: 25 pcs.

REF A.4021L CE / A.4021LST (sterile) CE Blue Scope Midi




Surgical/Examination Proctoscope trunk-end with led lighting built-in
 Tube length: 70mm
 Tube diameter: 21mm
 Blu Scope midi has similar design characteristics to code A.4023L/A.4023LST. Allows excellent control of the anal canal walls.
Quantity per box: 25 pcs.
 Also available in the version without adhesive
REF A.4021.1L.

REF A.4022L CE / A.4022LST (sterile) CE Blue Scope




Surgical/Examination Proctoscope trunk-end with led lighting built-in
 Tube length: 79mm
 Tube diameter: 23mm
 Blu Scope has similar design characteristics to code A.4023L/A.4023LST. Allows excellent control of the anal canal walls.
Quantity per box: 25 pcs.
 Also available in the version without adhesive
REF A.4022.1L.

REF A.4023L CE / A.4023LST (sterile) CE



Surgical/Examination Proctoscope flute beaked with led lighting built-in
 Tube length: 55/88mm
 Tube diameter: 23mm
 Lute beak to perform biopsies, elastic ligatures and haemorrhoid sclerosis, incision of intra-anal haematomas, excision of hypertrophic anal papillae, sutures with haemostatic stitiche following polypectomy, I.R. coagulation etc.
Quantity per box: 25 pcs.
 Also available in the version without adhesive
REF A.4023.1L.

REF A.4024L CE / A.4024LST (sterile) CE



Surgical/Examination Proctoscope trunk-end with led lighting built-in
 Tube length: 65mm
 Tube diameter: 23mm
 They have similar design characteristics to code A.4023L/A.4023LST. Allows excellent control of the anal canal walls.
Quantity per box: 25 pcs.
 Also available in the version without adhesive
REF A.4024.1L.