

A.4018

# Self·light®

sapi med  
DIAGNOSTIC  
LINE

*Disposable medical proctoscopes.*

A.4019

A.4021

A.4022

A.4023

A.4024



CE


The **Self·Light®** disposable proctoscope utilises a convenient battery-operated pen light, enabling examinations and minor therapeutic procedures to be performed in a non-specialist out-patients' department, on the ward or at home. **Self·Light®** can also be used with a plug-in light source produced by Sapi Med, if desired. **Self·Light®** proctoscopes are made of transparent plastic material and supplied in individually sealed bags. A free pen-light is included in each box of scopes.

**AVAILABLE IN STERILE VERSION AS WELL.**

## TECHNICAL FEATURES

- Clear, transparent plastic proctoscopes
- Flash-free
- Guaranteed atraumatic
- Compatible with all Sapi Med light sources or pen-light
- Conforms fully with essential requirements of European Council Directive 93/42/EU.
- Latex Free
- Non Sterile/Sterile


REF A.4018 CE / A.4018ST (sterile) CE



### Proctoscope diagnostic/paediatric flute beaked

Disposable, self-lighting  
 Tube length: 56mm  
 Tube diameter: 18mm  
 For clear and easy rectal examination.  
**Quantity per box:** 75 pcs. + 1 pen-light

REF A.4019 CE / A.4019ST (sterile) CE




### Proctoscope diagnostic/paediatric trunk-end

Disposable, self-lighting  
 Tube length: 54mm  
 Tube diameter: 18mm  
 For clear and easy rectal examination.  
**Quantity per box:** 75 pcs. + 1 pen-light

REF A.4021 CE / A.4021ST (sterile) CE

Blue Scope Midi




### Diagnostic/operative Proctoscope trunk-end

Disposable, self-lighting.  
 Tube length: 70mm  
 Tube diameter: 21mm  
**Box quantity:** 50 pcs.

REF A.4022 CE / A.4022ST (sterile) CE


Blue Scope



### Surgical/Examination Proctoscope trunk-end

Disposable, self-lighting.  
 Tube length: 79mm  
 Tube diameter: 23mm  
 It has an increase of about 1,4cm over the A.4024, so offering the practitioner improved access higher up the rectum. The introducer is designed for easy and atraumatic introduction to the patient.  
 A.4022 has similar design characteristics to code A.4023.  
 It allows excellent control of the anal canal walls.  
**Quantity per box:** 50 pcs. + 1 pen-light


REF A.4023 CE / A.4023ST (sterile) CE



### Surgical/Examination Proctoscope flute beaked

Disposable, self-lighting.  
 Tube length: 88mm maximum  
 Tube diameter: 23mm  
 Flute beak to perform biopsies, elastic ligatures and haemorrhoid sclerosis, incision of intra-anal haematomas, excision of hypertrophic anal papillae, sutures with haemostatic stitches following polypectomy, I.R. coagulation etc.  
**Quantity per box:** 50 pcs. + 1 pen-light

REF A.4024 CE / A.4024ST (sterile) CE



### Surgical/Examination Proctoscope trunk-end

Disposable, self-lighting.  
 Tube length: 65mm  
 Tube diameter: 23mm  
 A.4024 has similar design characteristics to code A.4023.  
 Allows excellent control of the anal canal walls.  
**Quantity per box:** 50 pcs. + 1 pen-light