

Schuco® Aspirators

Reliability you can count on.



S130A

S330A

S430A

Performance to meet your customers' toughest needs.

Schuco is widely known and respected in the hospital and homecare markets for reliability and durability. Now the entire family of Schuco aspirators has evolved to meet the changes in standards and provide just the right aspirator for your suction requirements.

- UL approved; meets IEC 60601-1
- Flow range higher than ever - over 40 lpm
- Every unit equipped with a hospital-grade plug
- Collection canisters with improved scratch and flake resistance
- High vacuum range for various suction needs - up to 560 mm Hg (22" Hg)
- Vibration-free gauge for easy reading and vacuum control
- Durable metal canister brackets accommodate a variety of collection containers
- 2-year warranty provides tangible proof of Schuco quality and durability

**Schuco**
byAllied

Schuco® Aspirators



S130A
14" L x 8" W x 11" H
(356 mm x 203 mm x 280 mm)
18 lbs. (8.16 kg)

S330A
15" L x 8" W x 10" H
381 mm x 203 mm x 254 mm
11.5 lbs. (5.21 kg)

S430A
14" L x 7" W x 10" H
(356 mm x 178 mm x 254 mm)
10.5 lbs. (4.76 kg)

SPECIFICATIONS

Vacuum Range:	Up to 560 mm Hg (22" Hg)
Flow Rate:	Greater than 40 lpm at regulator inlet
Sound Level:	Less than 58 dBA
Operating Temperature:	32°F-104°F (0°C-40°C)
UL Approval:	UL approved; meets IEC 60601-1
Warranty:	2-year limited warranty

ORDERING INFORMATION

S130A	Aspirator, 800 cc disposable canister, 1 tubing with filter, metal base, 115/60
S130A-GLASS	Aspirator, 750 cc reusable collection bottle, 1 tubing with filter, metal base, 115/60
S131A	Aspirator, 800 cc disposable canister, 1 tubing with filter, deluxe stand, 115/60
S330A	Aspirator, 800 cc disposable canister, 1 tubing with filter, molded base, 115/60
S430A	Aspirator, 800 cc disposable canister, 1 tubing with filter, 115/60, 4 legs (no base)

ACCESSORIES

S1160BA-RPL	Disposable Canister, 800 cc (5/case)
S1160BA-CS	Disposable Canister, 800 cc (70/case)
S1160-RPL	Disposable Canister, 1200 cc (4/case)
S1160A-RPL	Disposable Canister, 2000 cc (4/case)
S1161R	Reusable Collection Bottle, glass, 750 cc
S610100	Suction Tubing Kit (1 each 72" & 15" blue-tipped PVC tubing, 1/4" ID)
01-90-3928	Disposable Hydrophobic/Bacteria Filter, 1/4" barb to 1/4" barb (3/pkg.)
S1125	Stand for Schuco portable aspirators



S1125



S1125 shown
with S130A



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.
Toll Free: 800/444-3940
FAX: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
FAX: 800/246-6201

International
Telephone: 314/268-1683
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and
subject to change without notice.

Form 23-00-0001
Rev. January 2014